



**SklarTite™ Sterilization Container System**

# SklarTite™

- **LIGHT-WEIGHT  
AIRCRAFT ALUMINUM**
- **MINIMIZE STERILIZATION  
COSTS WITH PATENTED  
LOW PROFILE DESIGN**
- **30 DAY INTRODUCTORY  
PRICING FOR THIS INNOVATIVE  
NEW PRODUCT FROM SKLAR**

See reverse for ordering information

[www.sklarcorp.com](http://www.sklarcorp.com)

# SklarTite™ Sterilization Container System

Maximum Sterile Integrity • Economical, Low Cost-per-Cycle

The SklarTite™ line of sealed containers provides a simple yet effective system for the sterilization of operating room instruments. Designed for strength and light in weight, SklarTite™ containers are formed from an aircraft grade anodized aluminum alloy. Other sizes available by special order.



## SKLARTITE™ SEALED CONTAINERS

Perforated Top - Solid Bottom\*

### Extra-Small

7" Wide x 11" Long

10-1217 4" Height

### Small

12" Wide x 12" Long

10-1220 4" Height

10-1221 5" Height

### Medium

12" Wide x 17" Long

10-1226 4" Height

10-1227 5" Height

10-1228 6" Height

### Large

12" Wide x 23" Long

10-1233 5" Height

10-1234 6" Height



## SKLARTITE™ BASKETS

### Extra-Small Basket

6-1/2" Wide x 10" Long

10-1212 3" Height

### Small Basket

11" Wide x 11" Long

10-1242 3" Height

10-1244 4" Height

### Medium Basket

11" Wide x 16" Long

10-1250 3" Height

10-1252 4" Height


10-1253 5" Height

### Large Basket

11" Wide x 22" Long

10-1260 4" Height

10-1262 5" Height



## REPLACEMENT FILTERS, DISPOSABLE


1,000 Per Case

10-1278 7-3/4" x 7-1/4"

10-1279 9" x 9"

10-1281 9" x 6"

10-1283 9-1/2" x 9-1/2"



## SILICONE MATS

10-1282 Extra Small 6" x 10"

10-1284 Small 11" x 11"

10-1286 Medium 11" x 16"

10-1288 Large 11" x 22"



## FOAM INSERTS, DISPOSABLE

75 Per Pack, Blue

10-1207 3/16" x 20" x 10-1/2"

10-1205 3/8" x 20" x 10-1/2"

\* Perforated Bottoms Are Also Available.

Contact Your Sklar Account Manager For Details.

\*\* Other Sizes Available By Special Order

## FEATURES AND BENEFITS

### High Quality Construction

- Clear hard-coated aluminum
- Fade and Scratch-resistant
- Light weight, but high integrity
- Space efficient design
- Stronger than other systems

### Innovative cam-lock latching mechanism

- Positive gasket seal.

### Advanced filter system

- Rivetless, puncture proof filter protection
- Interchangeable filter plates
- Easy to use filters
- Field replaceable gaskets

### Unique Safety Enhancing Design

- Low edge, split shell design offers aseptic presentation protection that helps prevent accidental contamination of sets
- Space saving design requires only a 1-inch clearance between basket and container. Competitive products require 2 or more inches.

### Large Assortment of Accessories for organizing and Protecting General and Specialty/Surgical Instruments

- Silicone Mats
- Pins
- Dividers
- Cradling devices and posts



**Sklar**  
INSTRUMENTS

Phone: 610.430.3200

Toll Free: 800.221.2166

889 S. Matlack St.  
West Chester, PA 19382  
Surgi@sklarcorp.com  
www.sklarcorp.com