

SklarLite™ Rigid Sterilization Container System

Lightweight SklarLite™ Rigid Sterilization Containers provide an organized, routine and secure approach to sterilizing surgical instruments.

Features	Benefits
Lightweight anodized aluminum alloy	Ease of lifting and transportation of procedural sets.
Stackable Containers	Maximize space in autoclave and reduce cycle runs. Maximize storage shelf space and organize instrument sets by procedure.
Stackable Basket	Provide ease of instrument transfer per procedure, maximize space inside container during sterilization and eliminate the removal/disposal of wrapped tray material.
Instrument security	Delicate and laparoscopic instruments are secured in baskets/mats, minimizing damage to the instruments. Tamper evident seals prevent accidental sterile field breach.
Color Coded	Easy identification of instrument sets by procedure.
Ease of Filter Exchange	With a press of a button, filter is easily removed and replaced.
Safety Lid	Protection of sterile contents during transport/storage and ease of stacking.
Resusable	Minimize labor and eliminate material cost associated with wrapped trays.
270°F/4min. with 30min. Dry Time	Standard autoclave time/temperature; reduce cycle runs.

Portfolio

Available in Half, Mid and Full-Size, with depths from 4-10", containers come in your choice of silver, blue, green, red, yellow or black. Standard models consist of a bottom and perforated lid. Safe models consist of a bottom, perforated lid, and solid safety lid. Accessories sold separately and include wire baskets, silicone mats, disposable sterilization filters, process indicator data cards, ID labels and tamper evident seals.



	Half-Size			Mid-Size			Full-Size		
MODEL	NON-PERFORATED BOTTOM	PERFORATED BOTTOM	FEATURES	NON-PERFORATED BOTTOM	PERFORATED BOTTOM	FEATURES	NON-PERFORATED BOTTOM	PERFORATED BOTTOM	FEATURES
Standard	10-4500	10-4512	11.25" x 11" x 4"	10-4300	10-4310	18.25" x 11" x 4"	10-4000	10-4012	23" x 11" x 4"
Standard	10-4548	10-4560	11.25" x 11" x 5"	10-4348	10-4360	18.25" x 11" x 5"	10-4048	10-4060	23" x 11" x 5"
Standard	10-4596	10-4608	11.25" x 11" x 6"	10-4394	10-4406	18.25" x 11" x 6"	10-4096	10-4108	23" x 11" x 6"
Standard	10-4644	10-4656	11.25" x 11" x 8"				10-4144	10-4156	23" x 11" x 8"
Standard	10-4692	10-4704	11.25" x 11" x 10"				10-4192	10-4204	23" x 11" x 10"
Safe	10-4524	10-4536	11.25" x 11" x 4"	10-4324	10-4336	18.25" x 11" x 4"	10-4024	10-4036	23" x 11" x 4"
Safe	10-4572	10-4584	11.25" x 11" x 5"	10-4370	10-4382	18.25" x 11" x 5"	10-4072	10-4084	23" x 11" x 5"
Safe	10-4620	10-4632	11.25" x 11" x 6"	10-4418	10-4430	18.25" x 11" x 6"	10-4120	10-4132	23" x 11" x 6"
Safe	10-4668	10-4680	11.25" x 11" x 8"				10-4168	10-4180	23" x 11" x 8"
Safe	10-4716	10-4728	11.25" x 11" x 10"				10-4216	10-4228	23" x 11" x 10"