



F E M A L E P A T I E N T C A R E

# PESSARIES



- Popular Patterns
- Convenient Sizes
- Silicone, Latex Free
- Effective Treatment

[www.sklarcorp.com](http://www.sklarcorp.com)



# Pessaries

Sklar offers a full line of soft, silicone Pessaries for effective vaginal support and relief of incontinence. Pessaries are uniquely designed for ease of insertion and removal. Properly fitted, Pessaries can be successfully used to treat a variety of gynecologic conditions and abnormalities.

- Provides Effective Vaginal Support
- Soft, Silicone, Latex-Free
- Easy Insertion & Removal
- Includes Detailed Instructions

Choose from the most popular patterns, available in a variety of shapes and sizes. Latex Free.



PESSARY	PROLAPSE		CYSTOCELE	RECTOCELE
	1-2 Degree	3-4 Degree		
Ring with Support	X		X	
Donut		X	X	X
Gelhorn with Drain		X	X	
Cube with Drain		X	X	X

## RING WITH SUPPORT

For Mild Uterine Prolapse or a Cystocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5291	#1 2"	1 Each
90-5392	#2 2¼"	1 Each
90-5393	#3 2½"	1 Each
90-5294	#4 2¾"	1 Each
90-5295	#5 3"	1 Each
90-5296	#6 3¼"	1 Each
90-5397	#7 3½"	1 Each

## DONUT

For Moderate Uterine Prolapse or a Mild Cystocele or Rectocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5271	#1 2¼"	1 Each
90-5272	#2 2½"	1 Each
90-5273	#3 2¾"	1 Each
90-5274	#4 3"	1 Each
90-5275	#5 3¼"	1 Each
90-5276	#6 3½"	1 Each

## CUBE WITH DRAIN

For Moderate to Severe Uterine Prolapse or a Mild Cystocele or Rectocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5301	#1 29mm	1 Each
90-5302	#2 33mm	1 Each
90-5303	#3 37mm	1 Each
90-5304	#4 41mm	1 Each
90-5305	#5 45mm	1 Each

## GELHORN WITH DRAIN

For a Moderate Uterine Prolapse or Mild Cystocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5282	#2 2"	1 Each
90-5283	#3 2¼"	1 Each
90-5284	#4 2½"	1 Each
90-5285	#5 2¾"	1 Each
90-5286	#6 3"	1 Each

**For Price Quotes Contact Your Sklar Account Manager Or Distributor Sales Representative.**

USA

**800.221.2166**

INTERNATIONAL

**610.430.3200**

FAX

**610.696.9007**



889 S. Matlack St. West Chester, PA 19382 [www.sklarcorp.com](http://www.sklarcorp.com) [surgi@sklarcorp.com](mailto:surgi@sklarcorp.com) [international@sklarcorp.com](mailto:international@sklarcorp.com)