



S K L A R T E C H 5 0 0 0

LAPAROSCOPIC SYSTEMS

- Three Piece Modular System
- Miniature Instruments
- **NEW!** Single Port
- **NEW!** Disposables

www.sklarcorp.com

SINGLE PORT

NEW! THE NEXT EVOLUTION IN LAPAROSCOPIC SURGERY

Introducing a new line of single port instruments. Single port is an innovative procedure in which the skilled surgeon operates almost exclusively through a single entry point, typically the patients navel.

Patients will benefit from:

- Decreased trauma and risk for herniation.
- Enhanced cosmetic benefits with minimal visible scarring.
- Accelerated recovery time.
- Reduced post-operative pain.



SKLARTECH 5000

HIGHEST QUALITY... RELIABLE DURABILITY... EASE OF USE

Sklar offers a modular, three piece system of laparoscopic instruments and accessories that are the CHOICE of surgeons worldwide. With a complete line of the most popular patterns and interchangeable handles and replacement inserts, clinicians are guaranteed to find an instrument that meets their needs.

- Comfortable, carbon fiber, ratchet or non-ratchet ring handles.
- Surgical insert/tip is easily separated from cannula for cleaning & reprocessing.
- 360° Rotation allows for precise jaw positioning.
- Locking collar provides tension and prevents accidental movement.
- Insulated shaft protects against electrosurgical transfer.
- Miniature instruments available for delicate procedures.



DISPOSABLES

NEW! ECONOMIC SOLUTION WITHOUT COMPROMISING PERFORMANCE


Sklar offers your choice of two lines of disposable laparoscopic instruments. Both styles feature the most popular jaw patterns, excellent handle to blade responsiveness, insulated shaft to prevent arcing and burns, and ratchet or non-ratchet handle designs. With disposal of instruments immediately following procedures, the risk of cross-contamination, as well as reprocessing time and cost, are greatly reduced.

X-LINE

100% DISPOSABLE

HX LINE

DISPOSABLE SHEATH/INSERT COMBINED WITH REUSABLE HANDLE



**For Price Quotes Contact your Sklar Account Manager
or your Distributor Sales Representative for Details**

USA SALES

service@sklarcorp.com

PHONE: 610.430.3200

800.221.2166

INTERNATIONAL SALES

international@sklarcorp.com

PHONE: 610.430.3200



1333 Lenape Road
West Chester, PA 19382
www.sklarcorp.com