



RIGID STERILIZATION CONTAINER SYSTEM

**SKLARLITE™**

**50%  
LIGHTER!**



VISIT OUR WEBSITE [WWW.SKLARCORP.COM](http://WWW.SKLARCORP.COM) FOR MORE INFORMATION

# SKLARLITE™ RIGID STERILIZATION CONTAINER SYSTEM

**The SklarLite Rigid Sterilization Container System is a reusable rigid container made from anodized aluminum alloy and designed to be strong and durable yet lightweight.**

## STANDARD MODEL

Consists of a body and a lid and used with Sklar's disposable filters and other compatible accessories. Choose from perforated or non-perforated bottoms.

PRODUCT NUMBER		FEATURES
10-4500	Half Size Container	11.25" x 11" x 4"
10-4548	Half Size Container	11.25" x 11" x 5"
10-4596	Half Size Container	11.25" x 11" x 6"
10-4644	Half Size Container	11.25" x 11" x 8"
10-4300	Mid Size Container	18.25" x 11" x 4"
10-4348	Mid Size Container	18.25" x 11" x 5"
10-4394	Mid Size Container	18.25" x 11" x 6"
10-4000	Full Size Container	23" x 11" x 4"
10-4048	Full Size Container	23" x 11" x 5"
10-4096	Full Size Container	23" x 11" x 6"
10-4144	Full Size Container	23" x 11" x 8"

## SAFE MODEL

Consists of identical body and lid, but features an additional cover over the perforation field that safeguards contents during long term or external transport and storage.

PRODUCT NUMBER		FEATURES
10-4524	Half Size Container	11.25" x 11" x 4"
10-4572	Half Size Container	11.25" x 11" x 5"
10-4620	Half Size Container	11.25" x 11" x 6"
10-4668	Half Size Container	11.25" x 11" x 8"
10-4324	Mid Size Container	18.25" x 11" x 4"
10-4370	Mid Size Container	18.25" x 11" x 5"
10-4418	Mid Size Container	18.25" x 11" x 6"
10-4024	Full Size Container	23" x 11" x 4"
10-4072	Full Size Container	23" x 11" x 5"
10-4120	Full Size Container	23" x 11" x 6"
10-4168	Full Size Container	23" x 11" x 8"



Containers and accessories for the SklarLite Sterilization Container System are available in a variety of sizes and colors. Contact your Sklar Account Manager for information on the full line of products and customize a package that's right for you. Available colors: Silver, Yellow, Red, Blue, Green, Black.

**For Price Quotes Contact your Sklar Account Manager or your Distributor Sales Representative for Details**

- Pre-vacuum Steam
- Stackable Baskets
- Stackable Containers
- Transportable
- 6 month Shelf Life
- Cycle: 270°F/4min.
- 30 Minute Dry Time
- Load Limit 25lbs.

## ACCESSORIES

PRODUCT NUMBER		FEATURES
10-4788	Half Size Basket	10" x 9.5" x 3"
10-4790	Half Size Basket	10" x 9.5" x 4"
10-4482	Mid Size Basket	16" x 10" x 3"
10-4484	Mid Size Basket	16" x 10" x 4"
10-4288	Full Size Basket	21.25" x 10" x 3"
10-4290	Full Size Basket	21.25" x 10" x 4"
10-4970	Silicone Mat	8.5" x 9"
10-4975	Silicone Mat	15" x 9"
10-4985	Silicone Mat	20.5" x 9"
10-4900	Process Indicator Cards	250 / Bx
10-4910	Disposable Filters	100 / Pk
10-4925	Security Seals	100 / Pk
10-4930	Labels	Silver
10-4935	Labels	Yellow
10-4940	Labels	Red
10-4945	Labels	Blue
10-4950	Labels	Green
10-4955	Labels	Black



### USA SALES

[service@sklarcorp.com](mailto:service@sklarcorp.com)

PHONE: 610.430.3200

800.221.2166

### INTERNATIONAL SALES

[international@sklarcorp.com](mailto:international@sklarcorp.com)

PHONE: 610.430.3200



1333 Lenape Road  
West Chester, PA 19382  
[www.sklarcorp.com](http://www.sklarcorp.com)