



STERILE DISPOSABLE SAFETY SCALPELS

SKLARSAFE™



- “Click and Lock” Shield System
- You Never Touch the Blade!
- One-Hand Operation

www.sklarcorp.com



Sklarsafe™ Disposable Safety Scalpel

SklarSafe Sterile Disposable Safety Scalpels are conveniently equipped with movable sheaths to offer the surgeon protection and control without compromising performance or quality. Expose blade to allow for accurate device use. Cover blade to improve safety to health care workers and patients.

FEATURING THE UNIQUE

"CLICK AND LOCK" Safety Shield


- Protective, transparent sheath can be opened, closed, and locked with a single hand
- Audible "click" when sheath is placed in locked position
- Sturdy, "Sure-Grip" handle

Disposable Safety Scalpels are made with durable ABS plastic handles and fitted with surgical steel blades. Complies with current safety regulations for hospitals and surgery centers.

- Individually sealed in sturdy film and Tyvek® peel-back pouch
- Securely packaged in boxes of 10
- Sterilized with Gamma Radiation



06-3094 #10



06-3095 #11



06-3096 #15



SklarSafe™ Safety Scalpels

PN	Description	Quantities
06-3094	#10 Sterile	Pack of 10
06-3095	#11 Sterile	Pack of 10
06-3096	#15 Sterile	Pack of 10

**For Price Quotes Contact your Sklar Account Manager
or your Distributor Sales Representative for Details**

USA SALES

service@sklarcorp.com

PHONE: 610.430.3200

800.221.2166

INTERNATIONAL SALES

international@sklarcorp.com

PHONE: 610.430.3200



1333 Lenape Road
West Chester, PA 19382
www.sklarcorp.com

99-0172 02.2023
MKF910.05 Rev.01