



AUTOCLAVE AND INSTRUMENT CLEANER & CONDITIONER

SKLAR-SHEEN™



Sklar Instruments

SKLAR-SHEEN™

FOAM SPRAY
AUTOCLAVE CLEANER
& CONDITIONER

Conveniently Removes Rust,
Scale Spots, Stains & Corrosion

- Long Lasting Resonating Foam
Remains on inside walls of
autoclave

- Odor Counteractant
Eliminates unpleasant odors

- Rinses Freely
Imparts residue-free surface

- Non-Aerosol
Environmentally friendly

- Multi-Functional
Cleans & conditions high quality
stainless steel instruments,
carts & trays

Item # 10-2736 Contents: 22 fl. oz. (650 ml)


VISIT OUR WEBSITE WWW.SKLARCORP.COM

SKLAR-SHEEN™ is an Efficacious Phosphoric Acid Cleaning Agent ideal for surgical instruments and autoclaves.

EFFECTIVE IN LESS THAN 30 MINUTES

Sklar Sheen Foam Spray is an ideal cleaner for those hard to reach places. Convenient foam spray readily clings to autoclave walls. A robust rust and stain remover quickly cuts through steam sterilizer scale spots, stains, and deposits. It is also exceedingly effective at reconditioning high quality stainless steel instruments, carts, and trays.

- Restores Instruments to "Like-New" condition
- Reconditions Stainless Steel surfaces
- Extends Surgical Instrument Life
- Dissolves Hard Water Deposits
- Leaves Surfaces Residue-Free
- Eliminates Unpleasant Odors
- Frees Up Frozen Box Locks
- Fast and Effective
- Multi-Functional
- Non Corrosive
- Easy-to-Use



| PRODUCT NUMBER | DESCRIPTION | QUANTITY |
|----------------|---------------------------|----------|
| 10-2736 | 22 oz. Foam Spray Bottle | 1 Each |
| 10-2738 | 22 oz. Foam Spray Bottles | 12/case |

Sklar carries a full line of non-abrasive instrument and general purpose cleaning brushes to aid in the removal of stains, rust and debris from your stainless steel instruments and equipment.

| PRODUCT NUMBER | DESCRIPTION | QUANTITY | PRODUCT NUMBER | DESCRIPTION | QUANTITY |
|----------------|--|----------|----------------|---|----------|
| 10-1651 | Instrument Cleaning Brush, Nylon | 3/pk | 10-1422 | Autoclave Cleaning Brush with Extended Handle | 3/pk |
| 10-1418 | Cleaning Brush, 3/4" x 3" Brush Head | 3/pk | 10-1424 | Autoclave Cleaning Brush Set | Each |
| 10-1420 | Cleaning Brush, 7" x 1 1/4" Brush Head | 3/pk | 10-1434 | All Purpose Brush Set | Each |

For Price Quotes Please Contact Your Sklar Account Manager or Your Distributor Sales Representative.

USA

800.221.2166

INTERNATIONAL

610.430.3200

FAX

610.696.9007

889 S. Matlack St. West Chester, PA 19382 www.sklarcorp.com surgi@sklarcorp.com international@sklarcorp.com

99-0184 01.2014

