



# Sklar Soak™

Bactericidal • Virucidal • Fungicidal • Tuberculocidal

# Sklar Soak

For Professional Use • Ready To Use

**ACTIVE INGREDIENTS:**  
n-Alkyl (80% C14, 20% C16, 1% C12, 8% C18) dimethyl benzyl ammonium  
n-Alkyl (80% C12, 22% C14) ethylmethyl ethylbenzyl ammonium chloride  
Isopropylalcohol

**OTHER INGREDIENTS:**

TOTAL ..... CHILDREN

KEEP OUT OF  
CAUTION

See back panel for additional precautions.

**Contents: 1 U.S. Gallon / 3.785 Liters**

EPA REG. NO. 70144-1-31118

EPA EST. 1120-B-1

- ACTIVE INGREDIENTS:**  
n-Butyl 100%, CH<sub>3</sub>, 30% CH<sub>3</sub>, 13%, 8%, CH<sub>3</sub> dimethyl benzyl amine  
n-Butyl 100%, CH<sub>3</sub>, 32%, CH<sub>3</sub> dimethyl ethylbenzyl ammonium chloride  
Isopropyl  
Isopropyl
- OTHER INGREDIENTS:**  
TOTAL
- KEEP OUT OF REACH OF CHILDREN**
- CAUTION**
- See back panel for additional precautionary statements

**www.sklarcorp.com**






# Surgical Instrument Soak

## Sklar Soak™

### A Non-Glutaraldehyde Disinfectant Surgical Instrument Soak

Immersing soiled instruments in Sklar Soak™ for two minutes substantially reduces the risk of exposure to many infectious organisms including H1N1, HIV-1, HSV, VRE, Rotaviruses, Hepatitis B&C, as well as TB, MRSA, and A-Baumannii. If full immersion is not possible then cover soiled instruments with SKLAR SPRAY-ZYME enzymatic foam to prevent harmful encrustation while also precleaning.



**SKLAR INSTRUMENT SOAK™** is always effective as a decontaminant/disinfectant and this efficacy is without the irritating and caustic impact of glutaraldehydes, phenols and chlorine based products.

**SKLAR SOAK™** requires no special protection when handling nor are there any restrictive disposal guidelines.

- PH balanced, designed specifically for stainless steel instruments
- Intermediate level disinfectant with broad spectrum of efficacy
- Built in Rust Inhibitor
- Non-caustic
- Ready to use

## Sklar Soak™ Liquid

### Surgical Instrument Soak

Item #: 10-1603 Case of 4 Gallons

For Price Quotes and Samples Contact your Sklar Account Manager  
or your Distributor Sales Representative for Details

#### USA SALES

[service@sklarcorp.com](mailto:service@sklarcorp.com)

PHONE: 610.430.3200

800.221.2166

#### INTERNATIONAL SALES

[international@sklarcorp.com](mailto:international@sklarcorp.com)

PHONE: 610.430.3200



1333 Lenape Road  
West Chester, PA 19382  
[www.sklarcorp.com](http://www.sklarcorp.com)