



Sklar[®]

S U R G I C A L I N S T R U M E N T S O A K

Sklar Soak[™]

Sklar Instruments

Bactericidal • Virucidal* • Fungicidal • Tuberculocidal**

Sklar Soak

For Professional Use • Ready To Use

ACTIVE INGREDIENTS
n-alkyl 18% C14, 20% C16, 26% C18, 32% C20, 38% C22, 44% C24, 50% C26, 56% C28, 62% C30, 68% C32, 74% C34, 80% C36, 86% C38, 92% C40, 98% C42

OTHER INGREDIENTS

TOTAL

KEEP OUT OF REACH OF CHILDREN

CAUTION

See back panel for additional precautionary statements

Contents: 1 U.S. Gallon / 3.785 Liters

Cal. No. 10-1603

EPA REG. NO. 70144-1-31118

EPA EST. 1130-R-1

- **Non-Glutaraldehyde**
- **PH Neutral Formula Designed for Safe Use on All Surgical Instruments & Stainless Steelware**

www.sklarcorp.com



Surgical Instrument Soak

Sklar Soak™

A Non-Glutaraldehyde Product! Surgical Instrument Soak

Immersing soiled instruments in sklar soak™ for two minutes substantially reduces the risk of exposure to many infectious organisms including H1N1, HIV-1, HSV, VRE, Rotaviruses, Hepatitis B&C, as well as TB, MRSA, and A-Baumannii. If full immersion is not possible then cover soiled instruments with SKLAR SPRAY-ZYME enzymatic foam to prevent harmful encrustation while also precleaning.



SKLAR INSTRUMENT SOAK™ is always effective as a decontaminant/disinfectant and this efficacy is without the irritating and caustic impact of glutaraldehydes, phenols and chlorine based products.

SKLAR SOAK™ requires no special protection when handling nor are there any restrictive disposal guidelines.

- Non-Caustic Formula allows for use as a manual cleaner or in ultrasonic units.
- Effective Against: H1N1, HIV-1, HSV, VRE, Rotaviruses, Hepatitis B&C, as well as TB, MRSA, and A-Baumannii.
- EPA Registered
- Meets all OSHA standards on blood borne pathogen destruction
- Ready to use

Sklar Soak™ Liquid

Surgical Instrument Soak

Item# 10-1603 Case of 4 Gallons

For price quotes or samples, contact your sklar account manager or your distributor sales representative

USA SALES

PHONE: 610.430.3200
800.221.2166

INTERNATIONAL SALES

PHONE: 610.430.3200
FAX: 610.430.3941



889 S. Matlack St.
West Chester, PA 19382
surgi@sklarcorp.com
www.sklarcorp.com