



SURGICAL INSTRUMENT PRE-CLEANER

SPRAY-ZYME™



**PH Neutral Formula
Designed for Safe
Use on All
Instruments &
Stainless Steelware**

VISIT OUR WEBSITE AT WWW.SKILARCORP.COM

Multi-Enzyme Based Instrument Presoak & Cleaner

Foam Spray Cleaner

The newest product from Sklar for expert care of your surgical instruments. This multi-enzymatic spray creates a heavy foam which lasts for hours, keeping soiled instruments moist to avoid encrustation while also pre-cleaning prior to processing.

- Prevents Organic Encrustation - enzyme foam stops blood and bodily fluids from drying onto surgical instruments.
- No More Spills - messy spills between the operating room and central cleaning/processing department are eliminated.
- Neutral pH - safe and effective for stainless steel, tungsten carbide, aluminum and brass instruments.
- Bacteriostatic - microbiological control agents inhibit bacterial growth in the ready-to-use foam solution.
- Long Lasting Foam - saves time, foam lasts many hours after application
- Multi-Enzymatic - contains a synergistic blend of lipase and protease enzymes.



SKLAR SPRAY-ZYME

Item# 10-2722

Case of 12, 22 fl. oz.
Refillable Spray Bottles

Item# 10-2724

Case of 4, 1 Gallon

**For Price Quotes and Samples Contact your Sklar Account Manager
or your Distributor Sales Representative for Details**

USA SALES

service@sklarcorp.com

PHONE: 610.430.3200

800.221.2166

INTERNATIONAL SALES

international@sklarcorp.com

PHONE: 610.430.3200



1333 Lenape Road
West Chester, PA 19382
www.sklarcorp.com