



Surgical Instrument Pre-Cleaner

Spray-Zyme™



- **PH NEUTRAL FORMULA -
DESIGNED FOR SAFE USE ON
ALL INSTRUMENTS &
STAINLESS STEELWARE**

- **SPECIAL PRICING AT
WWW.SKILARCORP.COM**

www.sklarcorp.com

SKLAR SPRAY-ZYME

MULTI-ENZYME BASED INSTRUMENT PRESOAK & CLEANER



Foam Spray Cleaner

The newest product from Sklar for expert care of your surgical instruments. This multi-enzymatic spray creates a heavy foam which lasts for hours, keeping soiled instruments moist to avoid encrustation while also pre-cleaning prior to processing.

- *Prevents Organic Encrustation - enzyme foam stops blood and bodily fluids from drying onto surgical instruments.*
- *No More Spills - messy spills between the operating room and central cleaning/processing department are eliminated.*
- *Neutral pH - safe and effective for stainless steel, tungsten carbide, aluminum and brass instruments.*
- *Bacteriostatic - microbiological control agents inhibit bacterial growth in the ready-to-use foam solution.*
- *Long Lasting Foam - saves time, foam lasts many hours after application*
- *Multi-Enzymatic - contains a synergistic blend of lipase and protease enzymes.*

SKLAR SPRAY-ZYME

**Item# 10-2722 Case of 12, 22 fl. oz.
Refillable Spray Bottles**

Item# 10-2724 Case of 4, 1 Gallon

FOR PRICE QUOTES OR SAMPLES, CONTACT YOUR SKLAR ACCOUNT
MANAGER OR YOUR DISTRIBUTOR SALES REPRESENTATIVE

USA SALES

Phone: 610.430.3200
800.221.2166
Fax: 610.696.9007

INTERNATIONAL SALES

Phone: 610.430.3200
Fax: 610.430.3941
international@sklarcorp.com