



Carbon & Stainless Steel Blades / Disposable Scalpels

Sterile Disposable Blades & Scalpels

www.sklarcorp.com



STERILE DISPOSABLE BLADES AND SCALPELS

DISPOSABLE CARBON STEEL - Sterile, Box of 100

06-3010 #10 06-3012 #12 06-3020 #20 06-3022 #22
06-3011 #11 06-3015 #15 06-3021 #21 06-3023 #23

DISPOSABLE STAINLESS STEEL - Sterile, Box of 100

06-3100 #10 06-3102 #12 06-3104 #20 06-3106 #22
06-3101 #11 06-3103 #15 06-3105 #21 06-3107 #23

Securely and safely packaged in heavy foil and packed in a dispenser box of 100 blades.

DISPOSABLE SCALPELS - Sterile, Pack of 10

06-3110 #10 06-3112 #12 06-3120 #20 06-3122 #22
06-3111 #11 06-3115 #15 06-3121 #21 06-3123 #23

Securely and safely packaged in a peel back pouch and packed in boxes of 10.

Surgical Blades & Disposable Scalpels - Sklar No.3 & No.4 fitting surgical blades are manufactured in England in accordance with BS 2982 and BS EN 27740 in ISO 9000 certified facilities.

Sterile blades and disposable scalpels are sterilized by Gamma radiation to a minimum dose of 25 kGy.

Boxes are marked with lot number and expiration dates for traceability. All Sterile Product boxes are bar coded.

Microbiological tests are performed on all Sklar packaging materials and various other materials used in the manufacturing process. Bioburden estimates are randomly carried out on the materials and finished products from the Packing Department.

No post sterilization tests are performed as the product release is parametric. Testing for pyrogenicity is not required.

All Sklar precision products are latex free.

SKLAR NO.3 FITMENT SURGICAL SCALPEL BLADES

Fits Handles 3, 3L, 3 Graduated, 5B, 7, 9, B3 & B3L

Sterile Carbon Steel Surgical Blade Specifications

6, 9, 10, 10A, 11, 11P, E/11, Sabre E/11,12, 12d, 13, 14, 15, 15T, 15A, 15C, Sabre D/15, 16 and 40.

STEEL GRADE :	E1.30%
CARBON -	1.20% - 1.30%
SILICON -	0.10% - 0.35%
MANGANESE -	0.20% - 0.45%
CHROMIUM -	0.10% - 0.40%
SULPHUR -	0.025% MAX.
PHOSPHORUS -	0.035% MAX.

SKLAR NO.4 FITMENT SURGICAL SCALPEL BLADES

Fits Handles 4, 4L, 4 Graduated & 6B

Sterile Stainless Steel Surgical Blade Specifications

20, 21, 22, 22a, 23, 24, 25, 25A, 26 and 27.

STEEL GRADE :	SF.100
CARBON -	0.60% - 0.70%
SILICON -	0.50% MAX.
MANGANESE -	0.50% - 1.00%
CHROMIUM -	12.00% - 13.50%
SULPHUR -	0.025% MAX.
PHOSPHORUS -	0.03% MAX.

SKLAR DISPOSABLE SCALPELS

THE LEADING EDGE IN DISPOSABLE SCALPELS

Sklar disposable scalpels are produced with polystyrene handles and fitted with surgical steel blades from the No.3 & No.4 fitment ranges.

Steel grade and blade material analysis match the specifications of the sterile stainless steel surgical blades shown above.

"... Sklar Instruments had the blades I needed ... our supplier experienced a large back-order problem. '...upon delivery, I found the product to be superior and the customer service to be efficient."

Connie Dotson,
Baptist Hospital of
East Tennessee

Sklar Instruments was able to fulfill my blade requirements when no other supplier had the capability. My customer service representative was helpful and our service was prompt. In addition, Sklar was able to direct me to some other products which my staff is very satisfied with."

Sonny Wilson,
McKenzie Medical Center

Contact your Account Manager For Details

610.430.3200
800.221.2166
Surgi@sklarcorp.com
www.sklarcorp.com

