



Sterile, Disposable, Single-Use Biopsy Punch

Sklar Tru-Punch



- **SEAMLESS CUTTING EDGE FOR A PRECISE TISSUE SAMPLE**
- **CAREFULLY HONED EDGE FOR A SMOOTH RAZOR SHARP INCISION**
- **SPECIAL PRICING AVAILABLE AT WWW.SKILARCORP.COM**

www.sklarcorp.com



The **Sklar Tru-Punch** has a “seamless” steel tip carefully honed for a razor sharp incision and precise sample. The ribbed and hollow handle permits tissue to be easily dislodged.



SKLAR TRU-PUNCH

DISPOSABLE DERMAL BIOPSY PUNCH

BOX/25 STERILE

96-1097	1.5mm
96-1140	2mm
96-1142	3mm
96-1144	3.5mm
96-1146	4mm
96-1148	5mm
96-1150	6mm
96-1152	8mm
96-1160	Assorted - 5 of each: 2, 3, 3.5, 4, & 5mm



SKLAR TRU-PUNCH

DISPOSABLE DERMAL BIOPSY PUNCH

BOX/50 STERILE

96-1098	1.5mm
96-1100	2mm
96-1105	3mm
96-1110	3.5mm
96-1115	4mm
96-1120	5mm
96-1125	6mm
96-1130	8mm

**FOR PRICE QUOTES OR SAMPLES, CONTACT YOUR SKLAR ACCOUNT
MANAGER OR YOUR DISTRIBUTOR SALES REPRESENTATIVE**

USA SALES

Phone: 610.430.3200
800.221.2166
Fax: 610.696.9007

INTERNATIONAL SALES

Phone: 610.430.3200
Fax: 610.430.3941
international@sklarcorp.com



889 S. Matlack St.
West Chester, PA 19382
Surgi@sklarcorp.com
www.sklarcorp.com