



READY-TO-USE, PRE-MOISTENED TOWELETTES

# DISINFECTANT SURFACE WIPES



**Sklar Instruments**  
**Sklar Disinfectant**

## Surface Wipes

HOSPITAL/MEDICAL/SURGICAL/VETERINARY EQUIPMENT  
SURFACE DISINFECTANT AND CLEANER  
Hard, Non-Porous Surface Disinfectant Cleaner and Non-Food Contact Sanitizer  
Bactericidal • Fungicidal • Virucidal\* • Tuberculocidal

ACTIVE INGREDIENTS:	
n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	0.154%
dimethyl benzyl ammonium chloride	0.154%
n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> )	21.000%
dimethyl ethylbenzyl ammonium chloride	78.692%
Isopropanol	100.000%
OTHER INGREDIENTS	
TOTAL	

**KEEP OUT OF REACH OF CHILDREN**  
**CAUTION**

See side panels for First Aid and additional Precautionary Statements.

**Net Contents: 100 Wipes** **Net Weight: 1 lb 7.8 oz. (676.6 g)**

EPA REG. NO. 70144-2-31118 EPA EST. 56952-WI-001, 93725-TN-001

**Sklar Instruments**  
**Sklar Disinfectant**

## Surface Wipes

HOSPITAL/MEDICAL/SURGICAL/VETERINARY EQUIPMENT  
SURFACE DISINFECTANT AND CLEANER  
Hard, Non-Porous Surface Disinfectant Cleaner and Non-Food Contact Sanitizer  
Bactericidal • Fungicidal • Virucidal\* • Tuberculocidal

ACTIVE INGREDIENTS:	
n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	0.154%
dimethyl benzyl ammonium chloride	0.154%
n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> )	21.000%
dimethyl ethylbenzyl ammonium chloride	78.692%
Isopropanol	100.000%
OTHER INGREDIENTS	
TOTAL	

**KEEP OUT OF REACH OF CHILDREN**  
**CAUTION**

See side panels for First Aid and additional Precautionary Statements.

**Net Contents: 100 Wipes** **Net Weight: 1 lb 7.8 oz. (676.6 g)**

EPA REG. NO. 70144-2-31118 EPA EST. 56952-WI-001, 93725-TN-001

**2 Minute Kill Time for  
SARS-CoV-2 (Covid-19 Virus),  
MRSA AND E.COLI  
on hard, non-porous surfaces.**

**WWW.SKILARCORP.COM**





# SKLAR DISINFECTANT SURFACE WIPES

**Bactericidal • Virucidal • Fungicidal • Tuberculocidal**  
Ready-To-Use / Contains Biodegradable Detergent / Non-Bleach Formula

**2 Minute Kill Time for SARS-CoV-2 (Covid-19 Virus)**  
on hard, non-porous surfaces.

## 2 Minute Kill Time!

**SAME EFFECTIVE FORMULA AS  
SKLAR DISINFECTANT**

**CONVENIENT  
DISPENSING  
CANISTER**

- Ready-to-use, non-toxic formula
- **Kills 99.9% SARS-Related Coronavirus (SARS-CoV-2)** on hard, non-porous surfaces
- Specially formulated for the control of cross contamination
- Effective against H1N1, HIV-1, HSV, VRE, Rotaviruses, Hepatitis B&C, as well as TB, MRSA, E.Coli and A-Baumannii.
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas, laboratories and dialysis
- Satisfies OSHA's Bloodborne Pathogen Standards
- Destroys odor-causing bacteria
- 6" x 10" non-woven, tear resistant towelettes
- Clean scent
- Made with deionized water
- EPA registered

### SKLAR DISINFECTANT™ SURFACE WIPES

**Item # 10-1616 Case of 6**  
**100 Towelettes per Canister**

**For Price Quotes Contact your Sklar Account Manager  
or your Distributor Sales Representative for Details**

#### USA SALES

**service@sklarcorp.com**

**PHONE: 610.430.3200**

**800.221.2166**

#### INTERNATIONAL SALES

**international@sklarcorp.com**

**PHONE: 610.430.3200**



**1333 Lenape Road  
West Chester, PA 19382  
www.sklarcorp.com**

99-0173 05-2022  
MKF910.31 REV.06